



HoMedics Non-Contact Infrared Body Thermometer



Model: TE-350-EU

Manufacturer: HOMEDICS

EAN: 5010777151909

Product Code: 402910

R1,119

The HoMedics Non-Contact Infrared Body Thermometer allows for a non-invasive and non-contact temperature reading. Simply align the thermometer with the middle of the forehead and press the scan button to take a temperature measurement.

Non-Contact Temperature Measurement

The HoMedics Non-Contact Infrared Body Thermometer provides a clean and safe way to take multiple temperatures. The infrared sensor allows for non-contact measurements to be taken from 1-5cm away, offering extra protection against cross-contamination and spreading infections.

Instant Reading

The HoMedics Non-Contact Infrared Body Thermometer provides instant readings with a one-button control. After taking the measurement, the results will display in less than 2 seconds.

Easy To Read Icons

The product information provided is for reference purposes only. Technical specifications, photos, descriptions and prices may change without prior notice. The manufacturer and seller are not responsible for possible data discrepancies or typographical errors. We recommend checking the current information before purchasing.

itworkup.co.za | sales@itworkup.co.za | +27 (11)541-9950

After the measurement is complete, the temperature will be displayed on the LCD screen and the backlight colour will change to red, green or yellow. With the colour coded results, easily understand your measurements. Green – Normal Temperature. Yellow – Slight fever warning. Red – High fever.

Non-invasive

Using a standard thermometer can sometimes be quite invasive, particularly for children who aren't aware of why they're having their temperature taken from their ear, armpit or mouth. The HoMedics Non-Contact Infrared Body Thermometer is the ideal solution if you want to take readings without bothering your child. Also if you're monitoring the temperature of your baby, you won't need to disrupt them during their nap.

Store Up To 50 Temperature Readings

Easily keep track of your temperature readings thanks to the stored memory within the thermometer. Store up to 50 temperature readings.

To view your temperature history, simply press the MEMO button on the front of the device.

Surface Mode

The HoMedics Non-Contact Infrared Body Thermometer can also be used to check the temperature of foods and liquids such as baby milk. Also great for measuring the temperature of bath water.